

Utility of ICE COLD-PCR coupled to Sanger Sequencing or Next Generation Sequencing for increased mutation detection sensitivity during patient monitoring from circulating free DNA

Katherine Richardson², Filip Janku¹, Rui Lin², Ben Legendre², Gerald S. Falchook¹, Aung Naing¹, Veronica R. Holley¹, Siqing Fu¹, David S. Hong¹, Sarina A. Piha-Paul¹, Jennifer J. Wheeler¹, Ralph G. Zinner¹, Vivek Subbiah¹, Apostolia M. Tsimberidou¹, Daniel D. Karp¹, Vanda M. Stepanek¹, Goran Cabrilo¹, Rajyalakshmi Luthra³, Funda Meric-Bernstam¹, Agop Y. Bedikian⁴, Bryan K. Kee⁵, Cathy Eng⁵, Michael J. Overman⁵, Kevin B. Kim⁴, Amy Kruempel², Jaclyn Pope², Courtney Cubrich², Sarah Cherubin², Grant Wu², Vlada Melnikova², Marcia Lewis², Razelle Kurzrock⁶

¹Department of Investigational Cancer Therapeutics (Phase I Clinical Trials Program), MD Anderson Cancer Center, Houston TX; ²Transgenomic, Omaha, NE; ³Molecular Diagnostic Laboratory, ⁴Department of Melanoma Medical Oncology, ⁵Department of Gastrointestinal Medical Oncology, MD Anderson Cancer Center, Houston, TX; ⁶Moore's Cancer Center, San Diego, CA, USA

Introduction and Background:

Blood-based mutation analysis from serum and plasma is becoming very important for molecular diagnostics where no tumor is available as well as in the pharmacodynamic monitoring of the patient during therapy. Longitudinal monitoring from these liquid biopsies for the early detection of activating mutation recurrence or the emergence of resistance mutations is integral for ensuring proper treatment and best overall outcome for the patient. For monitoring purposes, a highly sensitive and rapid throughput technique is required.

ICE COLD-PCR (ICP) technology is capable of high sensitivity detection for both point mutations and indels through unbiased enrichment of relevant gene regions. Mutation detection in circulating cell-free DNA (cfDNA) is becoming increasingly important for monitoring treatment and resistance responses to various cancer drugs. Since the amount of cfDNA in plasma is limited, a methodology, such as ICP, that provides increased sensitivity for detecting these informative mutations is needed.

Next-Generation Sequencing (NGS) allows high throughput analysis of cancer mutations using target resequencing of a set/panel of genes, thus a broad mutation signature of tumors can be determined. However, the level of detection is ~2-4% for mutations unless a higher depth of coverage is used; this comes at the price of reducing the number of samples that can be analyzed on a chip. ICP enrichment of mutations prior to NGS mimics an increased "depth of coverage" without reducing the throughput per chip.

In this study, we show correlation data between solid tumor (standard analysis) and matched plasma using ICP coupled to Sanger sequencing. In addition, we show the advantages of using ICP as a pre-enrichment step for NGS.

Longitudinal Assessment of cfDNA BRAF & KRAS Mutations

Patient #	Disease	Tissue Mutation	Treatment prior to baseline cfDNA	Baseline cfDNA Mutations	Follow up cfDNA Mutation	Treatment Response*
5	Erdheim-Chester	BRAF V600E	No	BRAF V600E	Wild-Type	SD
6	Papillary Thyroid Cancer	BRAF V600E	No	BRAF V600E	Wild-Type	SD
2	Melanoma	BRAF V600E	No	BRAF V600E	Wild-Type	SD
3	Melanoma	BRAF V600E	No	BRAF V600E	Wild-Type	SD
4	Melanoma	BRAF V600E	No	BRAF V600E	Wild-Type	SD
1	Melanoma	BRAF V600E	Yes	Wild-Type	BRAF V600E	PD
7	Appendiceal Carcinoma	KRAS G13D	Yes	Wild-Type	KRAS G13D	PD
8	Colorectal Carcinoma	KRAS G12D	Yes	Wild-Type	KRAS G12D	PD


*SD = Stable disease; PD = Progressive disease

Concordance Study (ICP + Sanger): Tumor vs cfDNA, BRAF & KRAS

TESTED (N=31)	BRAF mutation CLIA	BRAF wild-type CLIA
BRAF mutation PLASMA	20	0
BRAF wild-type PLASMA	10	1
Observed agreements	21 (68%)	


TESTED (N=29)	KRAS mutation CLIA	KRAS wild-type CLIA
KRAS mutation PLASMA	22	2
KRAS wild-type PLASMA	4	1
Observed agreements	23 (79%)	

ICP-NGS (ICP prior to NGS analysis): Workflow



High Sensitivity Molecular Landscaping for Multiple Patient Samples


ICP-NGS: Proof-of-concept study




* More deletion reads than wild-type were filtered out by the variant caller.

Starting Percentage of Mutations in PIK3CA, KRAS and EGFR

Mutations enriched and detected individually using ICP-Sanger sequencing



NGS Analysis of cfDNA with and without ICE COLD-PCR Enrichment



Conclusions:

- ✓ Longitudinal Assessment Section:
 - ✓ ICP enrichment of mutations present in cfDNA provides:
 - ✓ A sensitive methodology for use in determining patient treatment options when tumor tissue is not available or difficult to obtain.
 - ✓ A monitoring capability to assess treatment effectiveness.
- ✓ Sanger Concordance Study Section:
 - ✓ ICP enrichment of mutations followed by Sanger sequencing show promising concordance with mutations detected in the tumor.
- ✓ ICP-NGS Section Workflow for multi-gene, high-throughput mutational analysis of mutations present in cfDNA:
 - ✓ ICP assays can be run individually and then combined and sequenced on NGS platforms.
 - ✓ ICP assays can be multiplexed and sequenced on NGS platforms.
- ✓ ICP-NGS Proof of Concept Section:
 - ✓ ICP prior to NGS analysis results in at least a 100-fold increase in the sensitivity of NGS alone.
 - ✓ Deletions are detected, but bioinformatics analyses for the deletions need refinement.
- ✓ NGS Analysis with and without ICP:
 - ✓ Enrichment of mutations by ICP prior to sequence analysis on NGS shows significant increases in mutation detection as compared to mutation detection using NGS alone.
- ✓ Next Steps:
 - ✓ Multiplex ICP for gene panels
 - ✓ Re-evaluate the concordance data using ICP to NGS.