


# BRAF and KRAS Mutation Testing in Plasma Cell-Free DNA with ICE COLD-PCR in Patients with Advanced Cancers (Selected Data from 2014 AACR Poster)

Filip Janku, Ben Legendre, Katherine Richardson, Gerald S. Falchook, Aung Naing, Veronica R. Holley, Siqing Fu, David S. Hong, Sarina A. Piha-Paul, Jennifer J. Wheler, Ralph G. Zinner, Vivek Subbiah, Apostolia M. Tsimberidou, Daniel D. Karp, Vanda M. Stepanek, Goran Cabrilo, Rajyalakshmi Luthra, Funda Meric-Bernstam, Agop Y. Bedikian, Bryan K. Kee, Cathy Eng, Michahel J. Overman, Kevin B. Kim, Amy Kruempel, Jaclyn Pope, Courtney Cubrich, Grant Wu, Marcia Lewis, Razelle Kurzrock

Investigational Cancer Therapeutics (Phase I Clinical Trials Program), The University of Texas MD Anderson Cancer Center, Houston, TX, USA; Transgenomic Inc., Omaha, NE; UCSD Moores Cancer Center, La Jolla, CA

## BACKGROUND

- Oncogenic mutations confer a survival and growth advantage to cancer cells
- Identifying oncogenic mutations in cancer can provide druggable targets for cancer therapies
- Plasma cell-free (cf) DNA in individuals with cancer offers an easily obtainable, low-risk, and inexpensive source of material for mutation analysis
- Longitudinal assessment of cfDNA can be used for monitoring of molecular changes throughout cancer therapy.



## METHODS


- Patients with advanced cancers, who were previously tested for *BRAF* V600 (42), *KRAS* G12/G13 (34), or both mutations (1) in the tumor samples (primary or metastatic) in a CLIA-certified laboratory during their clinical care were prospectively enrolled
- DNA from plasma (3-4ml) from patients with advanced cancers who progressed on systemic therapy were tested for *BRAF* V600 and *KRAS* G12 and G13 mutations using the ICE COLD-PCR platform
- ICE COLD-PCR, "Improved and Complete Enrichment COamplification at Lower Denaturation" selectively amplifies mutant DNA by exploiting differences in denaturation temperatures between mutant DNA duplexes and normal "wild-type" DNA duplexes
- *KRAS* Exon 2 and *BRAF* Exon 15 ICE COLD-PCR was performed on plasma samples
- Amplicons were analyzed by Sanger sequencing methods and results were compared to the mutation status of the archival primary or metastatic tumor tissue as determined in a CLIA-certified lab

## CONCLUSIONS


- ICE COLD-PCR detection of actionable mutations in *BRAF* and *KRAS* in cfDNA from plasma of patients with advanced cancers is feasible
- Longitudinal assessment of cfDNA mutations can demonstrate changes in mutation status during therapy, which seem to be in agreement with clinical course

## RESULTS: Longitudinal Assessment of cfDNA Mutations


**Patient 45: Metastatic Melanoma with *BRAF* V600E Mutation**



**Patient 22: Metastatic Melanoma with *BRAF* V600E Mutation**



**Patient 14: Metastatic Appendiceal Carcinoma with *BRAF* V600E Mutation**



**Patient 18: Metastatic Sigmoid Cancer with *KRAS* G13D Mutation**

