

CERTIFICATE OF ANALYSIS

Product: HemeScreen™ Cytopenia 3 Sample Plate

Product Number: HS-3P-CTP

Lot Number: 01310233-0922

Expiration Date: 2023-03

Precipio's HemeScreen™ Panel offerings are *in vitro* diagnostic assays that screen for clinically actionable somatic mutations in targets of interest related to chronic cytopenia. This reagent set is specifically designed to screen for mutations using High Resolution Melt Analysis in ASXL1, RUNX1, DNMT3A, SF3B1, and WT1 genes that are of clinical relevance.

Test	Specification	Pass/Fail
Functionality Tests	Amplification of controls.	Pass
	Amplification of genomic DNA.	Pass
	NTC free of target-specific amplification.	Pass
	Target sequence confirmation.	Pass
	Successful allelic discrimination between Wild-Type and Mutant Control.	Pass
Appearance	Desired wells contain HemeScreen™ reagents, a clear, colorless liquid with no particulate matter.	Pass
	Plates/Vials are sealed properly and show no sign of leakage.	Pass
	Plates/Vials inspected for any cracking after packaging.	Pass

Precipio, Inc.
 4 Science Park
 New Haven, CT USA
 Phone: 1-203-787-7888
www.precipiodx.com



Ayman Mohamed

Quality
 Issue Date: 9/25/2022